

JANSEN SELFRETAINING RETRACTOR

sharp; 3X3 teeth

10.0cm (4")




SKU: 070-630-10

Categories: [Open Surgery](#), [Retractors](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV030
KLS Martin	15-720-10
AUTHORS	JANSEN
Medicon	24.13.10
SPECIFICATIONS	sharp , toothed



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Selfretaining Retractors
Aesculap	BV030
KLS Martin	15-720-10

AUTHORS

[JANSEN](#)

Medicon

24.13.10

SPECIFICATIONS

[sharp, toothed](#)