

## JEWELER´S BIPOLAR COAGULATION FORCEPS

straight; blunt; Ø 0.6mm

11.0cm (4.4"); flat plug

**SKU:** 022-220-60

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.12.06
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.12.06
<b>SPECIFICATIONS</b>	<a href="#">blunt</a> , <a href="#">straight</a>