

JEWELER´S BIPOLAR COAGULATION FORCEPS

straight; blunt; Ø 0.6mm

11.0cm (4.4"); flat plug

SKU: 022-220-60

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.12.06
SPECIFICATIONS	blunt , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.12.06
SPECIFICATIONS	blunt , straight