

JEWELER´S BIPOLAR COAGULATION FORCEPS

angled; pointed; Ø 0.25mm

11.0cm (4.4"); flat plug

SKU: 022-221-25

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.12.07
SPECIFICATIONS	angled , sharp



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.12.07
SPECIFICATIONS	angled , sharp