

## JEWELER´S BIPOLAR COAGULATION FORCEPS

angled; pointed; Ø 0.25mm

11.0cm (4.4"); flat plug

**SKU:** 022-221-25

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>                           |
| <b>TYPE</b>           | <a href="#">Bipolar Coagulation Forceps</a>    |
| <b>Medicon</b>        | 88.12.07                                       |
| <b>SPECIFICATIONS</b> | <a href="#">angled</a> , <a href="#">sharp</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |  |
|-----------------------|--|
| <b>UNIT</b>           | <a href="#">pcs.</a>                           |
| <b>TYPE</b>           | <a href="#">Bipolar Coagulation Forceps</a>    |
| <b>Medicon</b>        | 88.12.07                                       |
| <b>SPECIFICATIONS</b> | <a href="#">angled</a> , <a href="#">sharp</a> |