

## JEWELER'S BIPOLAR COAGULATION FORCEPS

straight; pointed; Ø 0.25mm

11.0cm (4.4"); flat plug

**SKU:** 022-220-25

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.12.05
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.12.05
<b>SPECIFICATIONS</b>	<a href="#">sharp</a> , <a href="#">straight</a>