

JEWELER'S BIPOLAR COAGULATION FORCEPS

straight; pointed; Ø 0.25mm

11.0cm (4.4"); flat plug

SKU: 022-220-25

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Bipolar Coagulation Forceps |
| Medicon | 88.12.05 |
| SPECIFICATIONS | sharp , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Bipolar Coagulation Forceps |
| Medicon | 88.12.05 |
| SPECIFICATIONS | sharp , straight |