

## JEWELER´SBIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 0.6mm

11.0cm (4.4"); flat plug

**SKU:** 022-221-60

**Categories:** [Forceps](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.12.08
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">blunt</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bipolar Coagulation Forceps</a>
<b>Medicon</b>	88.12.08
<b>SPECIFICATIONS</b>	<a href="#">angled</a> , <a href="#">blunt</a>