

JEWELER´SBIPOLAR COAGULATION FORCEPS

angled; blunt; Ø 0.6mm

11.0cm (4.4"); flat plug

SKU: 022-221-60

Categories: [Forceps](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.12.08
SPECIFICATIONS	angled , blunt



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bipolar Coagulation Forceps
Medicon	88.12.08
SPECIFICATIONS	angled , blunt