

JOHNS-HOPKINS BULLDOG CLAMP

straight; length: 63mm


SKU: 026-018-06

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	13-119-06
AUTHORS	HOPKINS, JOHNS
Medicon	15.15.63
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
KLS Martin	13-119-06
AUTHORS	HOPKINS, JOHNS
Medicon	15.15.63
SPECIFICATIONS	straight