

JOHNS-HOPKINS BULLDOG CLAMP

straight; length: 75mm


SKU: 026-018-07

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	BH016
KLS Martin	13-119-07
AUTHORS	HOPKINS , JOHNS
Medicon	15.15.76
SPECIFICATIONS	straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Bulldog Clamps
Aesculap	BH016
KLS Martin	13-119-07
AUTHORS	HOPKINS , JOHNS
Medicon	15.15.76

SPECIFICATIONS

[straight](#)