

## JOHNS-HOPKINS BULLDOG CLAMP

straight; length: 75mm


**SKU:** 026-018-07

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	BH016
<b>KLS Martin</b>	13-119-07
<b>AUTHORS</b>	<a href="#">HOPKINS</a> , <a href="#">JOHNS</a>
<b>Medicon</b>	15.15.76
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	BH016
<b>KLS Martin</b>	13-119-07
<b>AUTHORS</b>	<a href="#">HOPKINS</a> , <a href="#">JOHNS</a>
<b>Medicon</b>	15.15.76

## SPECIFICATIONS

[straight](#)