

JOHNS-HOPKINS BULLDOG CLAMP

curved; length: 38mm

SKU: 026-019-04

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap BH011

AUTHORS [HOPKINS, JOHNS](#)

Medicon 15.15.38

SPECIFICATIONS [curved](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap BH011

AUTHORS [HOPKINS, JOHNS](#)

Medicon 15.15.38

SPECIFICATIONS [curved](#)