

JOHNS-HOPKINS BULLDOG CLAMP

curved; length: 63mm

SKU: 026-019-06

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Bulldog Clamps |
| Aesculap | BH014 |
| KLS Martin | 13-119-06 |
| AUTHORS | HOPKINS , JOHNS |
| Medicon | 15.15.63 |
| SPECIFICATIONS | curved |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|---|
| UNIT | pcs. |
| TYPE | Bulldog Clamps |
| Aesculap | BH014 |
| KLS Martin | 13-119-06 |
| AUTHORS | HOPKINS , JOHNS |
| Medicon | 15.15.63 |

SPECIFICATIONS

[curved](#)