

**JOHNS-HOPKINS BULLDOG CLAMP**

curved; length: 75mm

**SKU:** 026-019-07

**Categories:** [Hemostats](#), [Open Surgery](#)



**ZUSÄTZLICHE INFORMATIONEN**

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	BH015
<b>KLS Martin</b>	13-119-07
<b>AUTHORS</b>	<a href="#">HOPKINS</a> , <a href="#">JOHNS</a>
<b>Medicon</b>	15.15.70
<b>SPECIFICATIONS</b>	<a href="#">curved</a>



**ZUSÄTZLICHE INFORMATIONEN**

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	BH015
<b>KLS Martin</b>	13-119-07
<b>AUTHORS</b>	<a href="#">HOPKINS</a> , <a href="#">JOHNS</a>
<b>Medicon</b>	15.15.70

## SPECIFICATIONS

[curved](#)