

JOHNS HOPKINS BULLDOG CLAMP

straight; length: 38mm

SKU: 026-018-04

Categories: [Hemostats](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap BH011

AUTHORS [HOPKINS](#), [JOHNS](#)

SPECIFICATIONS [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Bulldog Clamps](#)

Aesculap BH011

AUTHORS [HOPKINS](#), [JOHNS](#)

SPECIFICATIONS [straight](#)