

## JOHNS HOPKINS BULLDOG CLAMP

straight; length: 38mm

**SKU:** 026-018-04

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	BH011
<b>AUTHORS</b>	<a href="#">HOPKINS, JOHNS</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Bulldog Clamps</a>
<b>Aesculap</b>	BH011
<b>AUTHORS</b>	<a href="#">HOPKINS, JOHNS</a>
<b>SPECIFICATIONS</b>	<a href="#">straight</a>