

JOHNS-HOPKINS HAEMOSTATIC CLAMP

curved; serrated

20.0cm (8")




SKU: 040-497-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Artery Forceps , Dissecting and Ligature Forceps |
| KLS Martin | 13-467-20 |
| AUTHORS | HOPKINS , JOHNS |
| Medicon | 17.07.45 |
| SPECIFICATIONS | curved , serrated |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|--|
| UNIT | pcs. |
| TYPE | Artery Forceps , Dissecting and Ligature Forceps |
| KLS Martin | 13-467-20 |
| AUTHORS | HOPKINS , JOHNS |
| Medicon | 17.07.45 |

SPECIFICATIONS

[curved](#), [serrated](#)