

JOHNS-HOPKINS HAEMOSTATIC CLAMP

curved; serrated

20.0cm (8")




SKU: 040-497-20

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
KLS Martin	13-467-20
AUTHORS	HOPKINS , JOHNS
Medicon	17.07.45
SPECIFICATIONS	curved , serrated



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Artery Forceps , Dissecting and Ligature Forceps
KLS Martin	13-467-20
AUTHORS	HOPKINS , JOHNS
Medicon	17.07.45

SPECIFICATIONS

[curved](#), [serrated](#)