

## JOHNSON NEEDLE HOLDER

curved; serrated; grooved jaw

27.0cm (10.8")




**SKU:** 050-193-27

**Categories:** [Open Surgery](#), [Suture](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>AUTHORS</b>	<a href="#">JOHNSON</a>
<b>Medicon</b>	10.43.87
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Needle Holders</a>
<b>AUTHORS</b>	<a href="#">JOHNSON</a>
<b>Medicon</b>	10.43.87
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">serrated</a>