

## JONES TOWEL FORCEPS

sharp

9.0cm (3.6")

**SKU:** 030-105-09

**Categories:** [Hemostats](#), [Open Surgery](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF407
<b>KLS Martin</b>	13-881-09
<b>AUTHORS</b>	<a href="#">JONES</a>
<b>Medicon</b>	17.50.08
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">sharp</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF407
<b>KLS Martin</b>	13-881-09

## AUTHORS

[JONES](#)

## Medicon

17.50.08

## SPECIFICATIONS

[curved](#), [sharp](#)