

JONES TOWEL FORCEPS

sharp; with tube holder

9.0cm (3.6")



SKU: 030-107-09

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
KLS Martin	13-883-09
AUTHORS	JONES
Medicon	17.60.08
SPECIFICATIONS	curved , sharp

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Towel Forceps
KLS Martin	13-883-09
AUTHORS	JONES
Medicon	17.60.08

SPECIFICATIONS

[curved](#), [sharp](#)