

JONES TOWEL FORCEPS

sharp; with tube holder

9.0cm (3.6")



SKU: 030-107-09

Categories: [Hemostats](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Towel Forceps |
| KLS Martin | 13-883-09 |
| AUTHORS | JONES |
| Medicon | 17.60.08 |
| SPECIFICATIONS | curved , sharp |

ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------|
| UNIT | pcs. |
| TYPE | Towel Forceps |
| KLS Martin | 13-883-09 |
| AUTHORS | JONES |
| Medicon | 17.60.08 |

SPECIFICATIONS

[curved](#), [sharp](#)