

## JONES TOWEL FORCEPS

sharp

5.0cm (2")



**SKU:** 030-105-05

**Categories:** [Hemostats](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF401
<b>KLS Martin</b>	13-881-05
<b>AUTHORS</b>	<a href="#">JONES</a>
<b>Medicon</b>	17.50.05
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">sharp</a>

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Towel Forceps</a>
<b>Aesculap</b>	BF401
<b>KLS Martin</b>	13-881-05

## AUTHORS

[JONES](#)

## Medicon

17.50.05

## SPECIFICATIONS

[curved](#), [sharp](#)