

SWING-SYSTEM ® KERRISON PUNCH

90° up cut; 3mm; with ejector

18.0cm (7.2"); 10mm open; CC



SKU: 249-051-03CC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| Medicon | 57.16.61 |
| SPECIFICATIONS | 90° up cut , with ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|-------------------------------------|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| Medicon | 57.16.61 |

SPECIFICATIONS

[90° up cut, with ejector](#)