

SWING-SYSTEM ® KERRISON PUNCH

90° up cut; 3mm; with ejector

18.0cm (7.2"); 10mm open; thin foot, CC




SKU: 249-051-03TCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	90° up cut , thin footplate , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	90° up cut , thin footplate , with ejector