

SWING-SYSTEM ® KERRISON PUNCH

90° up cut; 3mm; with ejector

18.0cm (7.2"); 10mm open; thin foot, CC



SKU: 249-051-03TCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--------------------------------------------------------------------------------------------|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 90° up cut , thin footplate , with ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--------------------------------------------------------------------------------------------|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 90° up cut , thin footplate , with ejector |