

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 3mm; with ejector

20.0cm (8"); 10mm open; CC



SKU: 249-072-03ACC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
Aesculap	FK914R
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
Aesculap	FK914R
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON

SPECIFICATIONS

[40° up cut, with ejector](#)