

## SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 4mm; with ejector

20.0cm (8"); 10mm open; CC



**SKU:** 249-072-04ACC

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>                                      |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>                       |
| <b>MATERIAL</b>       | <a href="#">Ceramic(CC)</a>                               |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>                                  |
| <b>Medicon</b>        | 999.34.19   |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">with ejector</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                 |                                     |
|-----------------|-------------------------------------|
| <b>UNIT</b>     | <a href="#">pcs.</a>                |
| <b>TYPE</b>     | <a href="#">Laminectomy Punches</a> |
| <b>MATERIAL</b> | <a href="#">Ceramic(CC)</a>         |
| <b>AUTHORS</b>  | <a href="#">KERRISON</a>            |
| <b>Medicon</b>  | 999.34.19                           |

## SPECIFICATIONS

[40° up cut, with ejector](#)