

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 4mm; without ejector

20.0cm (8"); 10mm open; CC



SKU: 249-072-04CC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
Aesculap	FK915M
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.71
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
Aesculap	FK915M
MATERIAL	Ceramic(CC)

AUTHORS

[KERRISON](#)

Medicon

57.17.71

SPECIFICATIONS

[40° up cut](#), [without ejector](#)