

## SWING-SYSTEM® KERRISON PUNCH

40° up cut; 2mm; with ejector; thin foot

23.0cm (9.2"); 10mm open




**SKU:** 249-073-02AT

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>AUTHORS</b>	<a href="#">KERRISON</a>
<b>SPECIFICATIONS</b>	<a href="#">40° up cut</a> , <a href="#">thin footplate</a> , <a href="#">with ejector</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Laminectomy Punches</a>
<b>AUTHORS</b>	<a href="#">KERRISON</a>
<b>SPECIFICATIONS</b>	<a href="#">40° up cut</a> , <a href="#">thin footplate</a> , <a href="#">with ejector</a>