

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 2mm; without ejector

23.0cm (9.2"); 10mm open; CC



SKU: 249-073-02CC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.32
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.32

SPECIFICATIONS

[40° up cut, without ejector](#)