

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 3mm; with ejector; thin foot

23.0cm (9.2"); 10mm open; CC




SKU: 249-073-03ATCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 40° up cut , thin footplate , with ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 40° up cut , thin footplate , with ejector |