

SWING-SYSTEM® KERRISON PUNCH

40° up cut; 3mm; without ejector

thin foot; 23.0cm (9.2"); 10mm open; CC



SKU: 249-073-03TCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| Medicon | 999.40.69 |
| SPECIFICATIONS | 40° up cut , thin footplate , without ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|-------------------------------------|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| Medicon | 999.40.69 |

SPECIFICATIONS

[40° up cut](#), [thin footplate](#), [without ejector](#)