

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 4mm; without ejector

thin foot; 23.0cm (9.2"); 10mm open; CC



SKU: 249-073-04TCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	999.40.70
SPECIFICATIONS	40° up cut , thin footplate , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	999.40.70

SPECIFICATIONS

[40° up cut](#), [thin footplate](#), [without ejector](#)