

## SWING-SYSTEM ® KERRISON PUNCH MODIFIZ.

40° up cut; 3mm; without ejector

thin foot; 20.0cm (8"); 15mm open; CC




**SKU:** 249-272-03TCC

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>  |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>   |
| <b>MATERIAL</b>       | <a href="#">Ceramic(CC)</a>   |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>  |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">thin footplate</a> , <a href="#">without ejector</a> |



### ZUSÄTZLICHE INFORMATIONEN

|                       |   |
|-----------------------|---|
| <b>UNIT</b>           | <a href="#">pcs.</a>  |
| <b>TYPE</b>           | <a href="#">Laminectomy Punches</a>   |
| <b>MATERIAL</b>       | <a href="#">Ceramic(CC)</a>   |
| <b>AUTHORS</b>        | <a href="#">KERRISON</a>  |
| <b>SPECIFICATIONS</b> | <a href="#">40° up cut</a> , <a href="#">thin footplate</a> , <a href="#">without ejector</a> |