

SWING-SYSTEM ® KERRISON PUNCH MODIFIZ.

40° up cut; 5mm; without ejector

thin foot; 20.0cm (8"); 15mm open; CC




SKU: 249-272-05TCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , thin footplate , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
SPECIFICATIONS	40° up cut , thin footplate , without ejector