

SWING-SYSTEM ® KERRISON PUNCH

90° up cut; 2mm; with ejector

18.0cm (7.2"); 10mm open; CC




SKU: 249-051-02CC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
Aesculap	FK949B
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.16.12
SPECIFICATIONS	90° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
Aesculap	FK949B
MATERIAL	Ceramic(CC)

AUTHORS

[KERRISON](#)

Medicon

57.16.12

SPECIFICATIONS

[90° up cut](#), [with ejector](#)