

SWING-SYSTEM ® KERRISON PUNCH

40° up cut; 1mm; without ejector

18.0cm (7.2"); 10mm open; CC




SKU: 249-071-01CC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.60
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.60

SPECIFICATIONS

[40° up cut, without ejector](#)