

SWING-SYSTEM ® KERRISON PUNCH

40° up; 2mm; without ejector

18.0cm (7.2"); 10mm open; CC



SKU: 249-071-02CC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.61
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.61

SPECIFICATIONS

[40° up cut, without ejector](#)