

SWING-SYSTEM ® KERRISON PUNCH


40° up cut; 2mm; without ejector

18.0cm (7.2"); 10mm open; thin foot, CC



SKU: 249-071-02TCC

Categories: [Neurosurgery](#), [Open Surgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.76
SPECIFICATIONS	40° up cut , thin footplate , without ejector

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
MATERIAL	Ceramic(CC)
AUTHORS	KERRISON
Medicon	57.17.76

SPECIFICATIONS

[40° up cut](#), [thin footplate](#), [without ejector](#)