

KERRISON PUNCH; 40° UP CUT

2mm, with ejector; thin footplate

23.0cm (9.2"); 10mm open; CC-coated



SKU: 240-073-02ATCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
AUTHORS	CLEANWAVE , KERRISON
SPECIFICATIONS	40° up cut , with ejector