

KERRISON PUNCH; 40° UP CUT

2mm; without ejector; thin footplate

18.0cm (7.2"); 10mm open; CC-coated



SKU: 240-071-02TCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 40° up cut , without ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 40° up cut , without ejector |