

KERRISON PUNCH; 40° UP CUT

5mm, with ejector; thin footplate

23.0cm (9.2"); 10mm open; CC-coated




SKU: 240-073-05ATCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| AUTHORS | CLEANWAVE , KERRISON |
| SPECIFICATIONS | 40° up cut , with ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|---|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| AUTHORS | CLEANWAVE , KERRISON |
| SPECIFICATIONS | 40° up cut , with ejector |