

KERRISON PUNCH; 40° UP CUT

3mm; without ejector

18.0cm (7.2"); 10mm open



SKU: 240-071-03

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
Aesculap	FF720
KLS Martin	24-993-03
AUTHORS	KERRISON
Medicon	57.08.43
SPECIFICATIONS	40° up cut , without ejector



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Laminectomy Punches
Aesculap	FF720
KLS Martin	24-993-03

AUTHORS

[KERRISON](#)

Medicon

57.08.43

SPECIFICATIONS

[40° up cut](#), [without ejector](#)