

KERRISON PUNCH; 40° UP CUT

3mm; without ejector

18.0cm (7.2"); 10mm open




SKU: 240-071-03

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| Aesculap | FF720 |
| KLS Martin | 24-993-03 |
| AUTHORS | KERRISON |
| Medicon | 57.08.43 |
| SPECIFICATIONS | 40° up cut , without ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-------------------|-------------------------------------|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| Aesculap | FF720 |
| KLS Martin | 24-993-03 |

AUTHORS

[KERRISON](#)

Medicon

57.08.43

SPECIFICATIONS

[40° up cut](#), [without ejector](#)