

KERRISON PUNCH; 40° UP CUT

2mm; with ejector

20.0cm (8"); 10mm open; CC-coated



SKU: 240-072-02ACC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 40° up cut , without ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 40° up cut , without ejector |