

KERRISON PUNCH; 90° UP CUT

2mm; with ejector; thin footplate

23.0cm (9.2"); 10mm open; CC-coated



SKU: 240-053-02TCC

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|------------------------------------------|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 90° up cut, with ejector |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|------------------------------------------|
| UNIT | pcs. |
| TYPE | Laminectomy Punches |
| MATERIAL | Ceramic(CC) |
| AUTHORS | KERRISON |
| SPECIFICATIONS | 90° up cut, with ejector |