

KIRSCHNER LIGATURE CONDUCTOR

strong curved

28.0cm (11.2")

SKU: 056-510-28

Categories: [Open Surgery](#), [Suture](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM924
KLS Martin	20-866-28
AUTHORS	KIRSCHNER
Medicon	13.30.28



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Ligature Conductors
Aesculap	BM924
KLS Martin	20-866-28
AUTHORS	KIRSCHNER

Medicon

13.30.28