

## KIRSCHNER LIGATURE CONDUCTOR

slightly curved

25.5cm (10")

**SKU:** 056-510-25

**Categories:** [Open Surgery](#), [Suture](#)



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Conductors</a>
<b>Aesculap</b>	BM925
<b>KLS Martin</b>	20-865-25
<b>AUTHORS</b>	<a href="#">KIRSCHNER</a>
<b>Medicon</b>	13.13.26
<b>SPECIFICATIONS</b>	<a href="#">curved</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Ligature Conductors</a>
<b>Aesculap</b>	BM925
<b>KLS Martin</b>	20-865-25

**AUTHORS**

[KIRSCHNER](#)

**Medicon**

13.13.26

**SPECIFICATIONS**

[curved](#)